

BLOOD TRANSFUSION ARTEFACT RECORD	
NUMBER	AA-028
ITEM / TITLE	A BOTTLE OF DRIED POOLED HUMAN PLASMA FOR TRANSFUSION
PHYSICAL DESCRIPTION	A 540 ml 'original style' Medical Research Council glass bottle of Dried Human Plasma; i.e. dried solids from 400 c.c. of pooled citrated plasma. The label is over-stamped: UNFILTERED. Some old cello tape remains around the metal cap.
OTHER INFORMATION	<p>The label on the bottle reads:</p> <p><i>Contains the dried solids (about 20 g. protein) from 400 c.c. of pooled citrated plasma, and may be given to patients of ANY BLOOD GROUP. To reconstitute add 400 c.c. sterile pyrogen-free distilled water. Intravenous saline or glucose solutions may also be used. Reconstituted plasma MUST BE DISCARDED if not used within 3 hours.</i></p> <p><i>Bottle No. K1307/2</i></p> <p><i>This number MUST be entered on the patient's notes and quoted in any enquiry.</i></p> <p>Note: The cello tape appears to have been used to help seal the bottle cap.</p>

